

How to use Fentanyl test strips

Fentanyl test strips can detect the presence of fentanyl in different types of drugs.

- ◆ Fentanyl test strips are not 100% accurate and drugs may still contain fentanyl. Use with caution.
- ◆ Fentanyl test strips indicate whether or not fentanyl may be present, but do not tell you how much fentanyl is present.
- ◆ Each test strip can only be used once. Discard after testing.



1. Prepare Supply

It is ideal to test all of your drugs. Test 5 mg of the substance (about the amount needed to cover half of Abe Lincoln's hair on a penny).

2. Mix with Water

Mix the powder into clean water (1 tsp) until it is completely dissolved. Fentanyl tends to clump together so be sure to mix thoroughly before testing.



3. Test Supply

Hold the solid blue end of the strip and dip the wavy line end of the strip in the water for 15 seconds.

4. Wait


Set the strip flat, away from anything absorbent. Wait 5 minutes.



How to Read Results

Positive Test:  Fentanyl Detected

Negative Test:  No Fentanyl Detected

Invalid Test:  Re-test with a New Strip